

## HYSTEOSALPINGOGRAPHY CATHETER



Reference:	Catheter size [F]	Fill liquid volume [ml]	Capacity of a balloon filled with air [cm <sup>3</sup> ]	Catheter length [cm]
KHSG5F	5	0.75	1.5	25
KHSG5F27				27
KHSG8F	8	2.5	3.5	25
KHSG8F27				27
KHSG12F	12	6.0	12	25

When ordering use Ref. Number, e.g.: KHSG5F.

Availability and specifications are subject to change. For more information, please contact your Balton sales representative.