

## DIALYSIS SET

double lumen



### Set components:

- double lumen catheter
- guidewire sleeve
- needle 18G, length 7 cm  
– for catheter size 6.5F  
needle 19G, length 7 cm
- guidewire
- dilator I
- dilator II

### Optional components



## Straight type

Reference: guidewire "J" type	Reference: guidewire Straight type	Catheter size [F]	Catheter length [cm]	Dilator I/II [F]/[F]	Guidewire size [inch]	Guidewire length [cm]	Version
ZDD6.5F10	ZDD6.5F10P	6.5	10	7/8	0.027"	60	Pediatric
ZDD6.5F12	ZDD6.5F12P		12				
ZDD7F10	ZDD7F10P	7	10	8/10	0.035"		
ZDD7F15	ZDD7F15P		15				
ZDD8F10	ZDD8F10P	8	10	10/12	0.035"		
ZDD8F12	ZDD8F12P		12				
ZDD8F15	ZDD8F15P		15				
ZDD9F10	ZDD9F10P	9	10	12/14	0.035"	For Adults	
ZDD9F15	ZDD9F15P		15				
ZDD11F15	ZDD11F15P	11	15				
ZDD11F18	ZDD11F18P		18				
ZDD11F20	ZDD11F20P		20				
ZDD11.5F15	ZDD11.5F15P	11.5	15				
ZDD11.5F20	ZDD11.5F20P		20				
ZDD12F15	ZDD12F15P	12	15				
ZDD12F18	ZDD12F18P		18				
ZDD12F20	ZDD12F20P		20				
ZDD13F15	ZDD13F15P	13	15				
ZDD13F18	ZDD13F18P		18				
ZDD13F20	ZDD13F20P		20				

When ordering use Ref. Number, e.g.: ZDD12F20.

### NOTE:

On demand, it is possible to select additional components.

While ordering, add the following endings to the reference index:

**S** – scalpel

**Y** – Y hub

**SR** – Guiding syringe 5 ml

**SL** – Syringe 5 ml

**KLLT** – TEGO connector

**KLLI** – Injection cap

**IV** – V needle

**N** – Nitinol guidewire

**H** – Hydrophilic coating

E.g.: ZDD12F20HPYKLLT.

Different dialysis catheter sizes and lengths available on demand. For detailed information, please contact the company's representative.

Availability and specifications are subject to change. For more information, please contact your Balton sales representative.

